

# CCSAP rabbit pAb

Cat No.:ES17759

For research use only


## Overview

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name             | CCSAP rabbit pAb                                                                                                                                                                                                                                                |
| Host species             | Rabbit                                                                                                                                                                                                                                                          |
| Applications             | WB                                                                                                                                                                                                                                                              |
| Species Cross-Reactivity | Human; Mouse                                                                                                                                                                                                                                                    |
| Recommended dilutions    | WB 1: 500-2000                                                                                                                                                                                                                                                  |
| Immunogen                | Synthesized peptide derived from human CCSAP AA range: 94-144                                                                                                                                                                                                   |
| Specificity              | This antibody detects endogenous levels of CCSAP at Human/Mouse                                                                                                                                                                                                 |
| Formulation              | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                         |
| Storage                  | Store at -20°C. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                              |
| Protein Name             | CCSAP                                                                                                                                                                                                                                                           |
| Gene Name                | CCSAP C1orf96 CSAP                                                                                                                                                                                                                                              |
| Cellular localization    | Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole . Cytoplasm, cytoskeleton, spindle . Cytoplasm, cytoskeleton . Cytoplasm, cytoskeleton, cilium basal body . Cytoplasm, cytoskeleton, cilium axoneme . Cell projection, axon . Cel |
| Purification             | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                           |
| Clonality                | Polyclonal                                                                                                                                                                                                                                                      |
| Concentration            | 1 mg/ml                                                                                                                                                                                                                                                         |
| Observed band            |                                                                                                                                                                                                                                                                 |
| Human Gene ID            | 126731                                                                                                                                                                                                                                                          |
| Human Swiss-Prot Number  | Q6IQ19                                                                                                                                                                                                                                                          |
| Alternative Names        |                                                                                                                                                                                                                                                                 |
| Background               |                                                                                                                                                                                                                                                                 |





**ELK Biotechnology**



Western blot analysis of lysates from SH-SY5Y cells,  
primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

[www.elkbiotech.com](http://www.elkbiotech.com)

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.