

BMP-8B rabbit pAb

Cat No.: ES1781

For research use only

Overview

Product Name BMP-8B rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human;Rat;Mouse;

Recommended dilutions Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not

yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human BMP8B. AA

range:261-310

Specificity BMP-8B Polyclonal Antibody detects endogenous

levels of BMP-8B protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Bone morphogenetic protein 8B

Gene Name BMP8B **Cellular localization** Secreted .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 45kD
Human Gene ID 656
Human Swiss-Prot Number P34820

Alternative Names BMP8B; BMP8; Bone morphogenetic protein 8B;

BMP-8; BMP-8B; Osteogenic protein 2; OP-2

Background This gene encodes a secreted ligand of the TGF-beta

(transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that

regulate gene expression. The encoded





+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com




preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer. The encoded protein stimulates thermogenesis in brown adipose tissue. Expression of this gene may be downregulated in pancreatic cancer. This gene may have arose from a gene duplication event and its gene duplicate is also present on chromosome 1. [provided by RefSeq, Jul 2016],





Western blot analysis of lysates from 293 and A549 cells, using BMP8B Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950

