



BCL9L (phospho-Ser915) rabbit pAb

Cat No.:ES18081

For research use only

Overview

Product Name	BCL9L (phospho-Ser915) rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized phospho peptide around human BCL9L (Ser915)
Specificity	This antibody detects endogenous levels of Human BCL9L (phospho-Ser915)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	BCL9L (Ser915)
Gene Name	BCL9L DLNB11
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	160kD
Human Gene ID	283149
Human Swiss-Prot Number	Q86UU0
Alternative Names	B-cell CLL/lymphoma 9-like protein (B-cell lymphoma 9-like protein) (BCL9-like protein) (Protein BCL9-2)

Background

