



C1qL2 rabbit pAb


Cat No.:ES1809

For research use only


Overview

Product Name	C1qL2 rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human C1QL2. AA range:231-280
Specificity	C1qL2 Polyclonal Antibody detects endogenous levels of C1qL2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Complement C1q-like protein 2
Gene Name	C1QL2
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	32-36kD
Human Gene ID	165257
Human Swiss-Prot Number	Q7Z5L3
Alternative Names	C1QL2; Complement C1q-like protein 2
Background	similarity:Contains 1 C1q domain.,similarity:Contains 1 collagen-like domain.,





Western Blot analysis of various cells using C1qL2 Polyclonal Antibody



Western blot analysis of lysates from HeLa and HUVEC cells, using C1QL2 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).

