

## ANR27 rabbit pAb

## Cat No.:ES18329

For research use only

## Overview


Product Name	ANR27 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse
<b>Recommended dilutions</b>	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human ANR27 AA
	range: 206-256
Specificity	This antibody detects endogenous levels of ANR27
	at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	ANR27
Gene Name	ANKRD27 PP12899
Cellular localization	Early endosome . Late endosome . Cytoplasmic
	vesicle membrane . Lysosome . Cell membrane .
	Melanosome . Colocalizes with VAMP7 in transport
	vesicles in the shaft of hippocampal neurons (By
	similarity)
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84079
Human Swiss-Prot Number	Q96NW4
Alternative Names	
Background	



ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Western blot analysis of lysates from HEK293 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C