

2B1E rabbit pAb

Cat No.: ES18539

For research use only

Overview

Product Name 2B1E rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human 2B1E AA

range: 194-244

Specificity This antibody detects endogenous levels of 2B1E at

Human

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name 2B1E Gene Name HLA-DRB1

Cellular localization

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band Human Gene ID

Human Swiss-Prot Number Q9GIY3


+86-27-59760950

Alternative Names

Background







Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night

