

1A36 rabbit pAb

Cat No.: ES18549

For research use only

Overview

Product Name 1A36 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human 1A36 AA

range: 150-200

Specificity This antibody detects endogenous levels of 1A36 at

Human

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name 1A36

Gene Name HLA-A HLAA

Cellular localization

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band Human Gene ID

Human Swiss-Prot Number P30455


Alternative Names

+86-27-59760950

Background







Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4° over night

