

1A1L2 rabbit pAb

Cat No.: ES18552

For research use only

Overview

Product Name 1A1L2 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human; Mouse Recommended dilutions WB 1:500-2000

Immunogen Synthesized peptide derived from human 1A1L2 AA

range: 453-503

Specificity This antibody detects endogenous levels of 1A1L2 at

Human/Mouse

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name 1A1L2 Gene Name ACCSL

Cellular localization

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 60kD
Human Gene ID 390110
Human Swiss-Prot Number Q4AC99


Alternative Names 1-aminocyclopropane-1-carboxylate synthase-like

protein 2 (ACC synthase-like protein 2)

Background







Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4° over night

