

# K1407 rabbit pAb

Cat No.:ES18864

For research use only

## Overview

|                          |   |
|--------------------------|---|
| Product Name             | K1407 rabbit pAb  |
| Host species             | Rabbit  |
| Applications             | WB  |
| Species Cross-Reactivity | Human;Rat;Mouse;  |
| Recommended dilutions    | WB 1:1000-2000  |
| Immunogen                | Synthesized peptide derived from human K1407 AA range: 720-770  |
| Specificity              | This antibody detects endogenous levels of human K1407  |
| Formulation              | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| Storage                  | Store at -20°C . Avoid repeated freeze-thaw cycles.   |
| Protein Name             | Coiled-coil domain-containing protein KIAA1407  |
| Gene Name                | KIAA1407  |
| Cellular localization    |   |
| Purification             | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality                | Polyclonal  |
| Concentration            | 1 mg/ml   |
| Observed band            |   |
| Human Gene ID            | 57577   |
| Human Swiss-Prot Number  | Q8NCU4  |
| Alternative Names        |   |
| Background               |   |

