

CNTD1 rabbit pAb

Cat No.:ES18867

For research use only

Overview

Product Name	CNTD1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CNTD1 AA range: 278-328
Specificity	This antibody detects endogenous levels of human CNTD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Cyclin N-terminal domain-containing protein 1
Gene Name	CNTD1 CNTD
Cellular localization	Nucleus . Cytoplasm . Chromosome . Shuttles between the nucleus and cytoplasm in a stage-specific manner of prophase I cells. Co-localized at crossover sites with PRR19. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	124817
Human Swiss-Prot Number	Q8N815
Alternative Names	
Background	

