

# EPMIP rabbit pAb

Cat No.:ES18870

For research use only

## Overview

Product Name	EPMIP rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human EPMIP AA range: 180-230
Specificity	This antibody detects endogenous levels of human EPMIP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	EPM2A-interacting protein 1 (Laforin-interacting protein)
Gene Name	EPM2AIP1 KIAA0766 My007
Cellular localization	Endoplasmic reticulum .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9852
Human Swiss-Prot Number	Q7L775
Alternative Names	
Background	

