



SRPRB rabbit pAb

Cat No.:ES18872

For research use only

Overview

Product Name	SRPRB rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SRPRB AA range: 174-224
Specificity	This antibody detects endogenous levels of human SRPRB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Signal recognition particle receptor subunit beta (SR-beta) (Protein APMCF1)
Gene Name	SRPRB PSEC0230
Cellular localization	Endoplasmic reticulum membrane ; Single-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	58477
Human Swiss-Prot Number	Q9Y5M8
Alternative Names	
Background	

