

TNPO3 rabbit pAb

Cat No.:ES18878

For research use only

Overview

Product Name	TNPO3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TNPO3 AA range: 601-651
Specificity	This antibody detects endogenous levels of human TNPO3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Transportin-3 (Importin-12) (Imp12) (Transportin-SR) (TRN-SR)
Gene Name	TNPO3 IPO12
Cellular localization	Nucleus envelope . Cytoplasm . Localizes to the nuclear envelope and annulate lamellae, which consists in stacks of endoplasmic reticulum membranes containing a high density of nuclear pores. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	23534
Human Swiss-Prot Number	Q9Y5L0
Alternative Names	
Background	

