



PRRT1 rabbit pAb

Cat No.:ES18882

For research use only

Overview

Product Name	PRRT1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PRRT1 AA range: 19-69
Specificity	This antibody detects endogenous levels of human PRRT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Proline-rich transmembrane protein 1 (Dispanin subfamily D member 1) (DSPD1)
Gene Name	PRRT1 C6orf31 NG5
Cellular localization	Cell membrane ; Single-pass type II membrane protein . Cell junction, synapse .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	80863
Human Swiss-Prot Number	Q99946
Alternative Names	
Background	

