



# XAF1 rabbit pAb

Cat No.:ES18885

For research use only

## Overview

<b>Product Name</b>	XAF1 rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB
<b>Species Cross-Reactivity</b>	Human;Mouse
<b>Recommended dilutions</b>	WB 1:1000-2000
<b>Immunogen</b>	Synthesized peptide derived from human XAF1 AA range: 47-97
<b>Specificity</b>	This antibody detects endogenous levels of human XAF1
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C. Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	XIAP-associated factor 1 (BIRC4-binding protein)
<b>Gene Name</b>	XAF1 BIRC4BP XIAPAF1
<b>Cellular localization</b>	Cytoplasm. Nucleus. Mitochondrion. Found in the cytoplasm and nucleus of placental syncytiotrophoblasts. Translocates to mitochondria upon TNF-alpha treatment.; [Isoform 1]: Nucleus.; [Isoform 5]: Nucleus.
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	
<b>Human Gene ID</b>	54739
<b>Human Swiss-Prot Number</b>	Q6GPH4
<b>Alternative Names</b>	
<b>Background</b>	

