



IGDC4 rabbit pAb

Cat No.:ES18886

For research use only

Overview

Product Name	IGDC4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human IGDC4 AA range: 56-106
Specificity	This antibody detects endogenous levels of human IGDC4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Immunoglobulin superfamily DCC subclass member 4 (Neighbor of punc e11) (Protein DDM36) (hDDM36)
Gene Name	IGDCC4 DDM36 KIAA1628 NOPE
Cellular localization	Cell membrane ; Single-pass type I membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	57722
Human Swiss-Prot Number	Q8TDY8
Alternative Names	
Background	

