



KRA15 rabbit pAb

Cat No.:ES18897

For research use only

Overview

Product Name	KRA15 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KRA15 AA range: 91-141
Specificity	This antibody detects endogenous levels of human KRA15
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Keratin-associated protein 1-5 (High sulfur keratin-associated protein 1.5) (Keratin-associated protein 1.5)
Gene Name	KRTAP1-5 KAP1.5 KRTAP1.5
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	83895
Human Swiss-Prot Number	Q9BYS1
Alternative Names	
Background	

