

RFOX3 rabbit pAb

Cat No.:ES18904

For research use only

Overview

Product Name	RFOX3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human RFOX3 AA range: 32-82
Specificity	This antibody detects endogenous levels of human RFOX3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	RNA binding protein fox-1 homolog 3 (Fox-1 homolog C)
Gene Name	RBFOX3
Cellular localization	Nucleus . Cytoplasm . Largely restricted to neuronal nuclei. However, significant cytoplasmic localization in neurons from brains from HIV-infected individuals with cognitive impairment. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	146713
Human Swiss-Prot Number	A6NFN3
Alternative Names	
Background	

