



# BRPF1 rabbit pAb

Cat No.:ES18906

For research use only

## Overview

Product Name	BRPF1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human BRPF1 AA range: 527-577
Specificity	This antibody detects endogenous levels of human BRPF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Peregrin (Bromodomain and PHD finger-containing protein 1) (Protein Br140)
Gene Name	BRPF1 BR140
Cellular localization	Nucleus . Chromosome . Cytoplasm . Localization to the nucleus depends on KAT6A, ING5 and MEAF6 (PubMed:18794358, PubMed:27939640). Localizes to transcription start sites (PubMed:24065767). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	7862
Human Swiss-Prot Number	P55201
Alternative Names	
Background	

