

# UB2J1 rabbit pAb

Cat No.:ES18907

For research use only

## Overview

Product Name	UB2J1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human UB2J1 AA range: 39-89
Specificity	This antibody detects endogenous levels of human UB2J1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Ubiquitin-conjugating enzyme E2 J1 (EC 6.3.2.19) (Non-canonical ubiquitin-conjugating enzyme 1) (NCUBE-1) (Yeast ubiquitin-conjugating enzyme UBC6 homolog E) (HsUBC6e)
Gene Name	UBE2J1 NCUBE1 CGI-76 HSPC153 HSPC205
Cellular localization	Endoplasmic reticulum membrane ; Single-pass type IV membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	51465
Human Swiss-Prot Number	Q9Y385
Alternative Names	
Background	

