



# K132L rabbit pAb

Cat No.:ES18913

For research use only

## Overview

Product Name	K132L rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human K132L AA range: 735-785
Specificity	This antibody detects endogenous levels of human K132L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	UPF0577 protein KIAA1324-like (Estrogen-induced gene 121-like protein) (hEIG121L)
Gene Name	KIAA1324L EIG121L
Cellular localization	Cell membrane ; Single-pass type I membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	222223
Human Swiss-Prot Number	A8MWY0
Alternative Names	
Background	

