

SEM3D rabbit pAb

Cat No.:ES18914

For research use only

Overview

Product Name	SEM3D rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SEM3D AA range: 442-492
Specificity	This antibody detects endogenous levels of human SEM3D
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	Semaphorin-3D
Gene Name	SEMA3D UNQ760/PRO1491
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	223117
Human Swiss-Prot Number	O95025
Alternative Names	
Background	

