

LYZL2 rabbit pAb

Cat No.:ES18919

For research use only

Overview

Product Name	LYZL2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LYZL2 AA range: 22-72
Specificity	This antibody detects endogenous levels of human LYZL2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Lysozyme-like protein 2 (Lysozyme-2) (EC 3.2.1.17)
Gene Name	LYZL2
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	119180
Human Swiss-Prot Number	Q7Z4W2
Alternative Names	
Background	

