



E4F1 rabbit pAb

Cat No.:ES18920

For research use only

Overview

Product Name	E4F1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human E4F1
Specificity	This antibody detects endogenous levels of human E4F1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Transcription factor E4F1 (EC 6.3.2.-) (E4F transcription factor 1) (Putative E3 ubiquitin-protein ligase E4F1) (Transcription factor E4F) (p120E4F) (p50E4F)
Gene Name	E4F1 E4F
Cellular localization	Nucleus, nucleoplasm. Cytoplasm. A small fraction is detected in the cytoplasm. Excluded from the nucleolus where it is targeted upon CDKN2A overexpression. Localizes to the mitotic spindle during embryogenesis. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	1877
Human Swiss-Prot Number	Q66K89
Alternative Names	
Background	

