

TCP1L rabbit pAb

Cat No.:ES18924

For research use only

Overview

Product Name	TCP1L rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TCP1L AA range: 68-118
Specificity	This antibody detects endogenous levels of human TCP1L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	T-complex protein 10A homolog 2 (T-complex protein 10A-2) (TCP10A-2) (TCP10-like)
Gene Name	TCP10L PRED77
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	140290
Human Swiss-Prot Number	Q8TDR4
Alternative Names	
Background	

