

DHX35 rabbit pAb

Cat No.:ES18928

For research use only

Overview

Product Name	DHX35 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DHX35 AA range: 618-668
Specificity	This antibody detects endogenous levels of human DHX35
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Probable ATP-dependent RNA helicase DHX35 (EC 3.6.4.13) (DEAH box protein 35)
Gene Name	DHX35 C20orf15 DDX35
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	60625
Human Swiss-Prot Number	Q9H5Z1
Alternative Names	
Background	

