



OR8U8 rabbit pAb

Cat No.:ES18929

For research use only

Overview

Product Name	OR8U8 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human OR8U8 AA range: 234-284
Specificity	This antibody detects endogenous levels of human OR8U8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Olfactory receptor 8U8
Gene Name	OR8U8
Cellular localization	Cell membrane; Multi-pass membrane protein.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	504189
Human Swiss-Prot Number	POC7N1
Alternative Names	
Background	

