

YIF1A rabbit pAb

Cat No.:ES18932

For research use only

Overview

Product Name	YIF1A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human YIF1A AA range: 101-151
Specificity	This antibody detects endogenous levels of human YIF1A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein YIF1A (54TMp) (YIP1-interacting factor homolog A)
Gene Name	YIF1A 54TM HYIF1P YIF1
Cellular localization	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Golgi apparatus membrane ; Multi-pass membrane protein . Endoplasmic reticulum-Golgi intermediate compartment membrane ; Multi-pass membrane protein . Cycles between the endoplasmic reticulum and the endoplasmic reticulum-Golgi intermediate compartment. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	10897
Human Swiss-Prot Number	O95070
Alternative Names	
Background	

