

METL5 rabbit pAb

Cat No.:ES18933

For research use only

Overview

Product Name	METL5 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human METL5 AA range: 59-109
Specificity	This antibody detects endogenous levels of human METL5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Methyltransferase-like protein 5 (EC 2.1.1.-)
Gene Name	METTL5 DC3 HSPC133
Cellular localization	Nucleus . Cell junction, synapse, presynapse . Cell junction, synapse, postsynapse .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	29081
Human Swiss-Prot Number	Q9NRN9
Alternative Names	
Background	

