

MB3L1 rabbit pAb

Cat No.:ES18935

For research use only

Overview

Product Name	MB3L1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MB3L1 AA range: 136-186
Specificity	This antibody detects endogenous levels of human MB3L1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Methyl-CpG-binding domain protein 3-like 1 (MBD3-like protein 1)
Gene Name	MBD3L1 MBD3L
Cellular localization	Nucleus . Nuclear, in large foci.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	85509
Human Swiss-Prot Number	Q8WWY6
Alternative Names	
Background	

