

SAMD9 rabbit pAb

Cat No.:ES18936

For research use only

Overview

Product Name	SAMD9 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SAMD9 AA range: 499-549
Specificity	This antibody detects endogenous levels of human SAMD9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Sterile alpha motif domain-containing protein 9 (SAM domain-containing protein 9)
Gene Name	SAMD9 C7orf5 DRIF1 KIAA2004 OEF1
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	54809
Human Swiss-Prot Number	Q5K651
Alternative Names	
Background	

