

# PCDG7 rabbit pAb

Cat No.:ES18937

For research use only

## Overview

Product Name	PCDG7 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PCDG7 AA range: 524-574
Specificity	This antibody detects endogenous levels of human PCDG7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protocadherin gamma-A7 (PCDH-gamma-A7)
Gene Name	PCDHGA7
Cellular localization	Cell membrane ; Single-pass type I membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	56108
Human Swiss-Prot Number	Q9Y5G6
Alternative Names	
Background	

