

## EX3L4 rabbit pAb

Cat No.:ES18939

For research use only

### Overview

Product Name	EX3L4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human EX3L4 AA range: 551-601
Specificity	This antibody detects endogenous levels of human EX3L4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Exocyst complex component 3-like protein 4
Gene Name	EXOC3L4 C14orf73
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	91828
Human Swiss-Prot Number	Q17RC7
Alternative Names	
Background	

