



REBL1 rabbit pAb

Cat No.:ES18940

For research use only

Overview

Product Name	REBL1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human REBL1 AA range: 32-82
Specificity	This antibody detects endogenous levels of human REBL1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	GTPase RhebL1 (Ras homolog enriched in brain like-1 c) (RhebL1c) (Ras homolog enriched in brain-like protein 1) (Rheb-like protein 1) (Rheb2)
Gene Name	RHEBL1
Cellular localization	Endomembrane system ; Lipid-anchor ; Cytoplasmic side . Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	121268
Human Swiss-Prot Number	Q8TAI7
Alternative Names	
Background	

