

TM10B rabbit pAb

Cat No.:ES18941

For research use only

Overview

Product Name	TM10B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TM10B AA range: 120-170
Specificity	This antibody detects endogenous levels of human TM10B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20 $^\circ\!\mathrm{C}$. Avoid repeated freeze-thaw cycles.
Protein Name	tRNA methyltransferase 10 homolog B (EC 2.1.1)
	(RNA (guanine-9-)-methyltransferase
	domain-containing protein 3)
Gene Name	TRMT10B RG9MTD3
Cellular localization	
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	158234
Human Swiss-Prot Number	Q6PF06
Alternative Names	
Background	



ELKbio@ELKbiotech.com

n www.elkb

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C