



CAND2 rabbit pAb

Cat No.:ES18946

For research use only

Overview

Product Name	CAND2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CAND2 AA range: 291-341
Specificity	This antibody detects endogenous levels of human CAND2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Cullin-associated NEDD8-dissociated protein 2 (Cullin-associated and neddylation-dissociated protein 2) (Epididymis tissue protein Li 169) (TBP-interacting protein of 120 kDa B) (TBP-interacting prote
Gene Name	CAND2 KIAA0667 TIP120B
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	23066
Human Swiss-Prot Number	O75155
Alternative Names	
Background	

