



LYZL4 rabbit pAb

Cat No.:ES18949

For research use only

Overview

Product Name	LYZL4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LYZL4 AA range: 88-138
Specificity	This antibody detects endogenous levels of human LYZL4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Lysozyme-like protein 4 (Lysozyme-4)
Gene Name	LYZL4 LYC4
Cellular localization	Secreted . Cytoplasmic vesicle, secretory vesicle, acrosome . Cell projection, cilium, flagellum . Found in the principal piece of sperm tail. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	131375
Human Swiss-Prot Number	Q96KX0
Alternative Names	
Background	

