

ARF3 rabbit pAb

Cat No.:ES18954

For research use only

Overview

Product Name	ARF3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ARF3 AA range: 66-116
Specificity	This antibody detects endogenous levels of human ARF3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	ADP-ribosylation factor 3
Gene Name	ARF3
Cellular localization	Golgi apparatus . Cytoplasm, perinuclear region .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	377
Human Swiss-Prot Number	P61204
Alternative Names	
Background	

