



EVX1 rabbit pAb

Cat No.:ES18955

For research use only

Overview

Product Name	EVX1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human EVX1 AA range: 312-362
Specificity	This antibody detects endogenous levels of human EVX1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Homeobox even-skipped homolog protein 1 (EVX-1)
Gene Name	EVX1
Cellular localization	Nucleus.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	2128
Human Swiss-Prot Number	P49640
Alternative Names	
Background	

