



SELO rabbit pAb

Cat No.:ES18959

For research use only

Overview

Product Name	SELO rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SELO AA range: 198-248
Specificity	This antibody detects endogenous levels of human SELO
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Selenoprotein O (SelO)
Gene Name	SELO
Cellular localization	Mitochondrion .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	83642
Human Swiss-Prot Number	Q9BVL4
Alternative Names	
Background	

