



SEHL2 rabbit pAb

Cat No.:ES18960

For research use only

Overview

| | |
|---------------------------------|---|
| Product Name | SEHL2 rabbit pAb |
| Host species | Rabbit |
| Applications | WB |
| Species Cross-Reactivity | Human;Rat;Mouse; |
| Recommended dilutions | WB 1:1000-2000 |
| Immunogen | Synthesized peptide derived from human SEHL2 AA range: 259-309 |
| Specificity | This antibody detects endogenous levels of human SEHL2 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C . Avoid repeated freeze-thaw cycles. |
| Protein Name | Serine hydrolase-like protein 2 (EC 3.1.-.-) |
| Gene Name | SERHL2 SERHL |
| Cellular localization | Cytoplasm, perinuclear region . Peroxisome . Concentrated in perinuclear vesicles. May be located in peroxisomes. . |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | |
| Human Gene ID | 253190 |
| Human Swiss-Prot Number | Q9H4I8 |
| Alternative Names | |
| Background | |

