

# HMN9 rabbit pAb

Cat No.:ES18962

For research use only

## Overview

Product Name	HMN9 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HMN9 AA range: 1-24
Specificity	This antibody detects endogenous levels of human HMN9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Humanin-like protein 9 (HN9) (MT-RNR2-like protein 9)
Gene Name	MTRNR2L9
Cellular localization	Secreted . Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	P0CJ76
Alternative Names	
Background	

