

HUNIN rabbit pAb

Cat No.:ES18963

For research use only

Overview

Product Name	HUNIN rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HUNIN AA range: 1-24
Specificity	This antibody detects endogenous levels of human HUNIN
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Putative humanin peptide
Gene Name	MT-RNR2
Cellular localization	Secreted . Cytoplasm . Cell projection, cilium, flagellum . Nucleus . Mitochondrion . Localizes to the sperm flagellum where it is highly concentrated in the midpiece (PubMed:20542501, PubMed:30920769). Detected in the cytoplasm and nucleus of spermatocytes and spermatids (PubMed:20542501). Also detected in sperm mitochondria (PubMed:20542501). In retinal pigment epithelium cells, detected in cytoplasm and mitochondria (PubMed:26990160). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q8IVG9
Alternative Names	





ELK Biotechnology

Background



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C