

ST65G rabbit pAb

Cat No.:ES18965

For research use only

Overview

Product Name	ST65G rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ST65G AA range: 191-241
Specificity	This antibody detects endogenous levels of human ST65G
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	STAGA complex 65 subunit gamma (Adenocarcinoma antigen ART1) (SPTF-associated factor 65 gamma) (STAF65gamma) (Suppressor of Ty 7-like)
Gene Name	SUPT7L KIAA0764
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9913
Human Swiss-Prot Number	O94864
Alternative Names	
Background	

