



TECT1 rabbit pAb

Cat No.:ES18968

For research use only

Overview

Product Name	TECT1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TECT1 AA range: 68-118
Specificity	This antibody detects endogenous levels of human TECT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Tectonic-1
Gene Name	TCTN1 TECT1 UNQ9369/PRO34160
Cellular localization	Cytoplasm, cytoskeleton, cilium basal body . Secreted . Despite the presence of a signal sequence, the full-length protein might not be secreted. Localizes at the transition zone, a region between the basal body and the ciliary axoneme. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	79600
Human Swiss-Prot Number	Q2MV58
Alternative Names	
Background	

