

PGRP2 rabbit pAb

Cat No.:ES18970

For research use only

Overview

Product Name	PGRP2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PGRP2 AA range: 90-140
Specificity	This antibody detects endogenous levels of human PGRP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	N-acetylmuramoyl-L-alanine amidase (EC 3.5.1.28) (Peptidoglycan recognition protein 2) (Peptidoglycan recognition protein long) (PGRP-L)
Gene Name	PGLYRP2 PGLYRPL PGRPL UNQ3103/PRO10102
Cellular localization	Secreted. Membrane.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	114770
Human Swiss-Prot Number	Q96PD5
Alternative Names	
Background	

