

## ZN483 rabbit pAb

Cat No.: ES18972

For research use only

## Overview

Product Name ZN483 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1:1000-2000

Immunogen Synthesized peptide derived from human ZN483 AA

range: 315-365

**Specificity** This antibody detects endogenous levels of human

ZN483

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Zinc finger protein 483 (Zinc finger protein with

KRAB and SCAN domains 16)

Gene Name ZNF483 KIAA1962 ZKSCAN16

Cellular localization Nucleus.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/ml

**Observed band** 

Human Gene ID 158399 Human Swiss-Prot Number Q8TF39

**Alternative Names** 

+86-27-59760950

**Background** 



23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.