

1A69 rabbit pAb

Cat No.: ES18976

For research use only

Overview

Product Name 1A69 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1:1000-2000

Immunogen Synthesized peptide derived from human 1A69 AA

range: 162-212

Specificity This antibody detects endogenous levels of human

1A69

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20℃. Avoid repeated freeze-thaw cycles.

Protein Name HLA class I histocompatibility antigen, A-69 alpha chain (Aw-69) (HLA class I histocompatibility antigen,

A-28 alpha chain) (MHC class I antigen A*69)

Gene Name HLA-A HLAA

Cellular localization

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/ml

Observed band Human Gene ID

Human Swiss-Prot Number P10316

Alternative Names

+86-27-59760950

Background

