



UK114 rabbit pAb

Cat No.:ES18981

For research use only

Overview

Product Name	UK114 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human UK114 AA range: 52-102
Specificity	This antibody detects endogenous levels of human UK114
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Ribonuclease UK114 (EC 3.1.-.-) (14.5 kDa translational inhibitor protein) (p14.5) (Heat-responsive protein 12) (UK114 antigen homolog)
Gene Name	HRSP12 PSP
Cellular localization	Cytoplasm . Nucleus . Peroxisome . Mitochondrion . Mostly cytoplasmic but, in less differentiated cells occasionally nuclear. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	10247
Human Swiss-Prot Number	P52758
Alternative Names	
Background	

