

# HDHD1 rabbit pAb

Cat No.:ES18986

For research use only

## Overview

Product Name	HDHD1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HDHD1 AA range: 89-139
Specificity	This antibody detects endogenous levels of human HDHD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Pseudouridine-5'-monophosphatase (5'-PsiMPase) (EC 3.1.3.n6) (Haloacid dehalogenase-like hydrolase domain-containing protein 1) (Haloacid dehalogenase-like hydrolase domain-containing protein 1A) (Pro
Gene Name	HDHD1 DXF68S1E FAM16AX GS1 HDHD1A
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	8226
Human Swiss-Prot Number	Q08623
Alternative Names	
Background	

